

Morning

LENS HANDLING

GUIDE



WASH HANDS WITH SOAP
AND WATER, THEN DRY

1.



INSERT 1-2 EYE
DROPS IN EACH EYE

2.



TAKE 3-4 HARD BLINKS.
FEEL THE LENS MOVE

3.



REMOVE THE **RIGHT**
CONTACT LENS

4.



CLEAN THE LENS IN THE
PALM OF YOUR HAND

5.



RINSE THE LENS WITH
DISINFECTION SOLUTION

6.



PLACE LENS IN THE CASE
SUBMERGE IN SOLUTION

7.



REPEAT SAME PROCESS
FOR **LEFT** CONTACT LENS

8.



CLOSE LENS CASE AND
STORE IN A SAFE PLACE

9.



Night

LENS HANDLING

GUIDE



GRAB YOUR CONTACT
LENS CASE




WASH HANDS WITH SOAP
AND WATER, THEN DRY



OPEN THE **RIGHT**
SIDE OF LENS CASE



INSPECT LENSES; NEVER
APPLY IF DIRTY OR DAMAGED



FILL BOWL OF LENS
WITH ARTIFICIAL TEARS



APPLY THE **RIGHT** CONTACT
TO YOUR RIGHT EYE



REPEAT STEP 4-6
FOR LEFT CONTACT LENS



EMPTY CONTENTS OF
LENS CASE INTO SINK



NEVER USE WATER
ON YOUR LENSES

